

AdVnt Biotechnologies, LLC

7898 E Acoma Drive, Suite 106 Scottsdale, AZ 85260

Web Catalog

Tularemia Threat Detection Kit 10/bx



Product URL <https://www.advntbiotechnologies.com/badd-tularemia-threat-detection-kit>

Short Description: AdVnt's BADD Tularemia Biowarfare Detection Test Kit- 10 Kits/Bx

Full Description: The BADD Tularemia Biowarfare Screening Assay is a Department of Homeland Security "Approved Product for Homeland Security" under the Support Anti-terrorism by Fostering Effective Technologies Act of 2002 (the SAFETY Act). This single-agent detection kit is designed to detect *Francisella tularensis* with a limit of detection (LOD) of 1.48×10^6 cfu/mL. Each kit contains everything needed to perform the test. No additional collection/dilution kits or readers are required. This item has a 2-year shelf life from date of manufacture and does not require refrigeration. Works in temperatures up to 120° Fahrenheit.

The BADD Tularemia detection kit includes the following:

- 1- pouched test
- 1- Instruction sheet
- 1- Buffer/Bottle/Dropper tube
- 1- Chain of Custody Label

Price: 270.00 USD

SKU: TUL-KIT-10

