



# AdVnt Biotechnologies, LLC

17200 W Bell Road, Surprise, AZ 85374

## Web Catalog

### PRO PAK COMBO "REFILL" 30/bx

Product URL <https://www.advntbiotechnologies.com/pro-pak-combo-refill-kit-pro-pak-combo-pro-strips-multi-agent-screening-system-10bx-plus-badd-anthrax-ricin-botulinum-plague-ypestis-seb-4eabx>



Short Description: Refill kit for the PRO PAK COMBO- 10 ProStrips (multi-agent) with 20 BADD (single-agent) Biological Agent Detection Devices.

Full Description: This is the refill kit for your PRO PAK Professional ballistic case. The PRO PAK COMBO REFILL is a combination of Advnt's Pro Strips Multi-agent and BADD single-agent BioWarfare Threat Detection screening assays. Both the "Pro Strips Multi-Agent and BADD Single-Agent Biowarfare Threat Screening Assay's are DHS "Approved Products for Homeland Security" under the Support Anti-terrorism by Fostering Effective Technologies Act of 2002 (the SAFETY Act). Pro Strips is the first multi-agent Biowarfare Threat Detection screening Kit that tests for five (5) biological targets [Anthrax, Ricin Toxin, Botulinum Toxin, Plague and SEB (Staphylococcal enterotoxin B)] with a single sample. BADD BioWarfare Threat Detection screening kits test for a single biological target with a single sample. The Pro Strip is in essence five BADD test strips in a five window, single sample HHA. The performance capabilities are the same for Pro Strips and BADD and are listed below: Threat Legal Limit of Detection Anthrax Bacillus anthracis;  $1.5 \times 10^4$  -  $8.3 \times 10^8$  /mL (strain dependent) Ricin Toxin Ricinus communis 10 ng/mL Botulinum Toxin Clostridium botulinum 33 ng/mL (BoNT A) 500 ng/mL (BoNT B) Plague Yersinia pestis  $1.0 \times 10^5$  CFU/mL SEB Staphylococcus aureus 10 ng/mL (Bacillus anthracis strains: Vollum, Ames, Sterne & New Hampshire) Each kit contains everything needed to perform the test. No additional Collection/Dilution kits or readers are required. This product has a 2 year (24 months) shelf life from date of manufacture and does not require special storage. Works in temperatures up to 120° Fahrenheit. Each PRO PAK REFILL includes the following: 20 - BADD single agent kits consisting of:

- 4 - Anthrax Threat Detection Kits
  - 4 - Ricin Toxin Threat Detection Kits
  - 4 - Botulinum Toxin Threat Detection Kits
  - 4 - Plague (Y. pestis) Threat Detection Kits
  - 4 - SEB Threat Detection Kits
  - 1 - recycleable shipping container
  - 1 - 8" x 11" laminated instruction card
- 10 - Pro Strips Rapid Screening Systems Multi-Threat Detection Devices
- Each Pro Strips & BADD Kit includes the following:
- 1- Test device
  - 1- Instruction sheet
  - 1- Buffer/swab/dropper tube
  - 1- Chain of Custody label

Price: 1410.00 USD

SKU: PS-5R30

