



AdVnt Biotechnologies, LLC

17200 W Bell Road, Surprise, AZ 85374

Web Catalog

Pro Strip Triple Toxin Environmental Threat Detection Kit 1/pouch



Product URL <https://www.advntbiotechnologies.com/pro-strip-triple-toxin-biowarfare-threat-detection-kit-1/pouch>

Short Description: Pro Strips PS-TOXIN Rapid Screening System Multi-toxin Threat Detection Kit 1/pouch

Full Description: Advnt's newly introduced TRIPLE TOXIN environmental threat detection screening kit that tests for three (3) biological targets with a single sample (Ricin, Botulinum and SEB). Toxin Causative Agent/Source LOD Botulinum Clostridium botulinum 33 ng/mL (BoNT A), 500 ng/mL (BoNT B) Ricin Ricinus Communis 10 ng/mL SEB Staphylococcus aureus 10 ng/mL Each kit contains everything needed to perform the test. No additional Collection/Dilution kits or readers are required. This product has a 2 year (24 months) shelf life from date of manufacture and does not require special storage. Made in the U.S.A./ Pro Strip Rapid Screening System's Commodity Jurisdiction Determination, by Department of State, is not subject to ITAR jurisdiction. The Department of Commerce advises the item is classified as EAR99. Each Pro Strips Triple-Toxin- Rapid Screening System Kit includes the following:

- 1- Multi-agent test device
- 1- Instruction sheet
- 1- Buffer/swab/dropper tube
- 1- Chain of Custody label

Price: 75.00 USD

SKU: PS-TOXIN

